




Erkopress Ci Motion





Erkopress Ci Motion: The Ci Motion is additionally equipped with an integrated powerful compressor with reserve compressed-air.

- Pressure Molding Technique
- Infrared Temperature Sensor
- No Preheating Time.
- Vertical Forming
- Touchscreen Operation
- Automated Process (Motorized)
















Occluform-3

The Erkoform-3D Motion: The unit is based on vacuum forming technique. This unit consists of an automated thermoforming process.

Automation removes the need to monitor the heating process and manually moving the foil holder from the heating element to the model to thermoform.

3D + :Thermoforming unit with inbuilt reserve vacuum and touch-panel program operation the same characteristics as the Erkoform-3dmotion but without automated thermoforming process.

Occluform-3: The Occluform-3 can be installed on the Erkoform 3D+ and 3D Motion unit. The opposing bite can be articulated with occluform on the thermoformed sheet to register the occlusal contacts.



